

GFP-Atg8 Cleavage Assay

Hitoshi Nakatogawa, Hayashi Yamamoto

METHOD

For yeast strain construction, culture, and induction of non-selective autophagy by starvation or rapamycin treatment, also see the protocols "Cultivation of Yeast Cells and Induction of Autophagy" and "Fluorescence Microscopy of GFP-Atg8 Transport into the Vacuole and Observation of Autophagic Bodies".

Sample preparation

1. Culture yeast cells expressing GFP-Atg8 to mid-log phase ($OD_{600} \approx 1.5$) at 30°C, and induce autophagy as described in the protocol "Cultivation of Yeast Cells and Induction of Autophagy".
2. Take 1 ml of the culture, and harvest cells by centrifugation at 15,000 *g* for 1 min at 4°C using a centrifuge (MX-105, TOMY). Store the cell pellets at -30°C until all samples are obtained.

NOTE GFP fragments generated by vacuolar degradation of GFP-Atg8 become detected in immunoblotting analysis after nitrogen starvation/rapamycin treatment for ~2 h under our experimental conditions. We usually take samples at 0, 2, and 4 h.

NOTE An *atg* mutant should be included as a control in order to judge whether GFP-Atg8 cleavage occurs in an autophagy-dependent manner under your experimental conditions.

NOTE Cell pellets can be stored at -30°C for a few weeks.

3. Suspend the cells in 100 μ l of ND solution. Stay on ice for 15 min.
4. Mix with 100 μ l of 20% TCA. Stay on ice for 15 min.
5. Centrifuge at 15,000 *g* for 5 min at 4°C.
6. Remove the supernatants, and suspend the pellets in 0.5 ml of cold acetone using a bath sonicator containing iced water.
7. Centrifuge at 15,000 *g* for 5 min at 4°C.
8. Remove the supernatants, and dry the pellets at room temperature for 10 min with the lids of the test tubes opened.
9. Dissolve the pellets in [the value of OD_{600} units \times 50] μ l of SDS sample buffer by mixing at 65°C for 10 min using a mixer (Thermomixer Comfort, Eppendorf),

and then boil for 3 min.

NOTE Some proteins, at least Ape1, can disappear in immunoblotting detection, if the pH of the sample buffer becomes acidic due to the carry-over of TCA. Thus, we recommend using SDS-PAGE sample buffer of the composition described below.

10. Centrifuge at 15,000 *g* for 1 min at room temperature, and use the supernatants for immunoblotting analysis.

SDS-PAGE and Immunoblotting

11. To prepare 10 ml of 12% acrylamide gel, mix 4 ml of 30% acrylamide/bisacrylamide (37.5:1), 2.5 ml of 1.5 M Tris-HCl (pH 8.8), 100 μ l of 10% (w/v) SDS, and 3.3 ml of H₂O.

12. Add 50 μ l of 20% (w/v) APS and 4 μ l of TEMED, and gently mix.

13. Cast the mixture into assembled glass plates (10 cm \times 16 cm \times 1 mm).

14. Gently overlay a few milliliters of 70% EtOH. Stay at room temperature for 15 min.

15. Remove 70% EtOH thoroughly.

16. Mix 2.5 ml of Stacking gel mix with 12.5 μ l of 20% (w/v) APS and 2.5 μ l of TEMED, overlay the mixture onto the separate gel, insert a comb, and stay at room temperature for 15 min.

17. Remove the comb and flush lanes with distilled water.

18. Load 10 μ l of the samples and molecular weight markers (Precision Plus Protein™ Dual Color Standards, Bio-Rad, 1610374), and perform electrophoresis at a constant current of 40 mA for ~90 min using Gel electrophoresis apparatus (BE-260/BE-250, Bio Craft, and PowerPac HC, BioRad), and Running buffer for SDS-PAGE (Nacalai Tesque) (until the dye front reaches the bottom of the gel).

19. Set filter papers (Blot Paper, Bio Craft), a PVDF membrane (Immobilon-P PVDF membranes, Millipore), and the gel on a semidry transfer cell (Trans-Blot Turbo System, Bio-Rad), and perform transfer at a constant voltage of 25 V for 25 min using Transfer buffer.

20. Shake the membrane in Blocking buffer at room temperature for 30 min using a seesaw shaker (BC-700, Bio Craft).

21. Exchange Blocking buffer and add anti-GFP antibody to a dilution of 1:2,000. Shake at 4°C overnight.
22. Shake the membrane in TBS-T for 5 min at room temperature. Repeat this wash step three times.
23. Shake the membrane in Blocking buffer containing peroxidase-conjugated rabbit anti-mouse IgG (1:5,000 dilution) at room temperature for 30 min.
24. Wash the membrane in TBS-T for 5 min at room temperature three times.
25. Detect the bands of GFP-Atg8 and GFP fragments on the membrane using chemiluminescent reagent for horseradish peroxidase (FEMTOGLOW^{Plus}, Michigan Diagnostics LLC) and a imager system (LAS4010 ImageQuant System, GE Healthcare). The bands of GFP-Atg8 and GFP fragments will appear near the molecular weight markers of 37 kD and 25 kD, respectively.

MATERIALS

REAGENTS

- NaOH (Wako, 198-13765)
- Dithiothreitol (DTT) (Nacalai Tesque, 14112-52)
- Trichloroacetic acid (TCA) (Nacalai Tesque, 34637-85)
- Acetone (Nacalai Tesque, 00310-53)
- Tris(hydroxymethyl)aminomethane (Tris) (Nacalai Tesque, 35434-21)
- HCl (Nacalai Tesque, 18321-05)
- Sodium dodecyl sulfate (SDS) (Nacalai Tesque, 31607-65)
- Glycerol (Nacalai Tesque, 17018-83)
- Bromophenol blue (Wako, 021-02911)
- 30% acrylamide and bisacrylamide (37.5:1) solution (Bio-Rad, 06144-05).
- Ammonium persulfate (APS) (Nacalai Tesque, 02627-34)
- N, N, N', N' -tetramethylethylenediamide (TEMED) (Nacalai Tesque, 33401-72)
- Glycine (Wako, 073-00737)
- Precision Plus Protein™ Dual Color Standards (Bio-Rad, 1610374)
- Running Buffer Solution (10x) for SDS-PAGE, Tris-Glycine (Nacalai Tesque, 30329-74)
- Methanol (Nacalai Tesque, 21915-93)
- Immobilon-P PVDF membranes (Millipore, IPVH00010).

- Semidry Gel Pad (filter paper) (Bio Craft, 7002)
- Grade 3MM Chr Blotting Paper (GE Healthcare, 3030-917)
- NaCl (Wako, 191-01665)
- KCl (Wako, 163-03545)
- Tween-20 (Nacalai Tesque, 28353-85)
- Meiji Step skim milk (Meiji)
- Anti-GFP monoclonal antibody (Clontech, JL-8, 632381)
- Peroxidase-conjugated rabbit anti-mouse IgG (Jackson ImmunoResearch, 315-035-003)
- FEMTOGLOW^{Plus} (Funakoshi/Michigan Diagnostics LLC, SHRPE21008)

EQUIPMENT

- Centrifuge MX-105 (TOMY)
- Thermomixer Comfort (Eppendorf)
- Gel electrophoresis apparatus BE-260, BE-250 (Bio Craft)
- PowerPac HC (BioRad)
- Trans-Blot Turbo System (Bio-Rad)
- Seesaw Shaker BC-700 (Bio Craft).
- LAS4010 ImageQuant System (GE Healthcare)

REAGENT PREPARATION

ND solution

- 0.2 M NaOH
- 0.1 M DTT

To prepare 1 ml of this reagent, mix 0.2 ml of 1 M NaOH, 0.1 ml of 1 M DTT (store at -20°C), and 0.7 ml of distilled water, and chill on ice, immediately before use.

SDS sample buffer

- 0.1 M Tris-HCl (pH 7.5)
- 2% (w/v) SDS
- 10% (v/v) glycerol
- 20 mM DTT

- bromophenol blue (trace amount)

To prepare 50 ml of this buffer, mix 5 ml of 1 M Tris–HCl (pH 7.5) with 10 ml of 10% (w/v) SDS, 5 ml of glycerol, 0.154 g of dithiothreitol, 30 ml of distilled water, and a trace amount of bromophenol blue. Store aliquots at -20°C .

Stacking gel mix

- 5% acrylamide and bisacrylamide (37.5:1) solution
- 125 mM Tris–HCl (pH 6.8)
- 0.1% (w/v) SDS

To prepare 300 ml of this buffer, mix 50 ml of 30% acrylamide and bisacrylamide solution (37.5:1) with 75 ml of 0.5 M Tris–HCl (pH 6.8), 3 ml of 10% (w/v) SDS and 169 ml of ultrapure distilled. Store at 4°C .

Transfer buffer

Lower Buffer (for top thick papers)

- 0.3 M Tris
- 10% (v/v) methanol

To prepare 1 L of this buffer, mix 200 ml of 1.5 M Tris, 100 ml of methanol, and 700 ml of ultrapure distilled. Store at room temperature.

Middle Buffer (for middle thin papers and PVDF membranes)

- 25 mM Tris
- 0.02% (w/v) SDS
- 10% (v/v) methanol

To prepare 1 L of this buffer, mix 16.8 ml of 1.5 M Tris, 100 ml of methanol, 2 ml of 10% (w/v) SDS, and 881.2 ml of distilled water. Store at room temperature.

Upper Buffer (for bottom thick papers)

- 25 mM Tris
- 0.02% (w/v) SDS
- 10% (v/v) methanol
- 40 mM 6-aminohexanoic acid

To prepare 1 L of this buffer, mix 16.8 ml of 1.5 M Tris, 100 ml of methanol, 2 ml of 10% (w/v) SDS, 5.2 g of 6-aminohexanoic acid, and 881.2 ml of distilled water. Store at room temperature.

Put papers and PVDF membranes in the following order (from the bottom to the top): a thick paper, a thin paper, a PVDF membrane, a gel, and a thick paper.

TBS-T (10×)

- 250 mM Tris
- 138 M NaCl
- 27 mM KCl
- 0.01% (v/v) Tween 20

To prepare 1 L of this 10× stock solution, mix 30.25 g of Tris, 80 g of NaCl, 2 g of KCl, and 10 ml of Tween-20 in 900 ml of distilled water. Add 12.5 ml of 1 M HCl to adjust the pH to 7.4, and top up to 1 L with distilled water. For a 1× solution, mix 100 ml of this stock solution with 900 ml of distilled water.

Blocking buffer

- 2% skim milk in TBS-T buffer

To prepare 50 ml of this buffer, dissolve 1 g of skim milk in TBS-T.